



Pill Buddie Setup Guide (short)

DISCLAIMER: This unit is *NOT* a medical device. It is a vending device designed to dispense pills, vitamins, supplements and even candy. If used for any type of medication products, this should be done by a qualified person. Please do not add any liquids to this device.

It is recommended to wear gloves or wash your hands prior to managing any medication. Always use the **Stylus** when touching the screen. NEVER use your fingers or any sharp devices such as a pen.

Step-by-Step User Instructions

Step 1: Choose a Location

1. Place the Pill Buddie on a level work surface, near a power outlet.

Step 2: Connect the Power

1. Plug the power adapter into the back of the device, and connect to a wall outlet.
2. Press the **PWR** (Power) button on the back of the device.

Step 3: Wi-Fi Setup

NOTE: The very first time the startup completes, or if you move to a new network, Pill Buddie will prompt you to setup your Wi-Fi.

1. Use the stylus to select the  **Wi-Fi** icon.
2. Pill Buddie will search for all local networks. If you don't see yours, tap **Scan** to refresh the list.
3. Tap the drop-down to choose your **Network**.
4. Tap the password field and enter your **Password**.
5. Tap the **Connect** button join the desired **Network**.

Step 4: Time Zone/DST, Sound and Vibration Setup

1. Use the stylus to select the  **Setup** icon.
2. Tap the drop-down menu to choose the **Time Zone** from your area.
3. To enable sound, enable the optional **Sound On** check box.
4. If your area follows **Daylight Savings Time (DST)**, enable the checkbox. If not, leave unchecked.
5. If enabled, the **Vibrator** option helps to prevent larger pills, vitamins or candy from getting stuck in the carousel.

Step 5: Loading Pills, Vitamins or Candy

1. Collect and sort all tablets prior to being added.
2. Use the tweezers to place each one into the correct slot in the carousel. Do not overfill!
3. Each slot contents depend on the time-of-day pills, vitamins or candy are added.

. For example, if you have the following events:

- 9:00 AM



- 8:00 PM

If you load the Pill Buddie **before** 9:00 AM or **after** 8:00 PM, then you can start loading with the morning medication in slot 1.

If you load the Pill Buddie **after** 9:00 AM and **before** 8:00 PM then you must load the Pill Buddie with the second medication in slot 1 since the next scheduled event is 8:00 PM, the evening dose.

NOTE: If pills, vitamins or candy are added at noon, the first dispensation will be the next PM time. Likewise, if pills, vitamins or candy are added after the last PM event, dispensing will be at the next AM event.

Step 6: Set Your Dispensing Schedule

1. Use the stylus to tap the  **Schedule** screen.
2. Use the stylus to select the Event to enable.
 - Select the hour, minute and AM/PM from the popup window.
 - Select Save

Step 7: Collecting Your Pills, Vitamins or Candy

1. Each time pills, vitamins or candy are dispensed, blue LEDs are displayed. 
 - If **Sound On** and **Vibrator** options are enabled, you will hear a sound, and the Pill Buddy vibrates.
2. Take your pills, vitamins or candy from the tray.
 - Removing the tray turns the blue LEDs red .
3. Return the tray to continue the dispensing schedule. The LED lights will go out until the next dispensation time.

Need Help? For support or questions, contact:

PillBuddie Healthy Tech

 info@PillBuddie.ca

 (902) 989-1807

 <https://pillbuddie.ca>